



BELOTERO® SOFT

Product Specifications

1-844-226-8277

info@medsupplysolutions.com

medsupplysolutions.com

BELOTERO®

MERZ AESTHETICS

SOFT LIDOCAINE

Caution

Consult instructions for use

Do not use if package is damaged

Single use product. Do not re-use

Open the blister by pulling the Tyvek lid following the arrow

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.

Sterile. Sterilized by irradiation. Only the needle itself is sterile, but not the outside of the needle packaging.

Temperature limit of storage: 2 °C – 25 °C

Batch code

Use-by date

CE mark in accordance with Directive 93/42/EEC relating to medical devices. This mark is followed by the notified body number.

Date of manufacture

Manufacturer

Manufacturer of the needles:

TSK Laboratory, Japan,
2-1-5 Hirayanagi-cho, Tochigi-Shi,
Tochigi-Ken, 328-0012 Japan;

The needles are CE marked  0123

Manufacturer of Belotero Soft Lidocaine:








ANTEIS SA
18 Chemin des Aulx
CH-1228 Plan-les-Ouates
Geneva, Switzerland

Australian Sponsor Name and Address:





Merz Australia Pty Ltd
Level 3, 244 Coward Street
Mascot, NSW, 2020
Australia

Date of the instructions for use: 2021-11-05

13. Instructions how to use the patient implant card
An implant card is provided with Belotero Soft Lidocaine. This implant card must be completed by the physician according to the below instructions and provided to the patient after injection.

Number	Symbols	Details
1		Please enter patient's name.
2		Please enter the date of implantation.
3-4		Please enter the name and address of the healthcare professional.
5		Please enter the number of injections.
6		Please enter the total volume injected.
7		Please enter the injection site(s)
8		Please stick here one of the product traceability label, you can keep the second one for your records.

Important information to be provided to the patient
Instruct the patient to keep the patient implant card with her/him and to present it to her/his healthcare professional in case of other appointments. Information about previous treatment must be presented to her/his healthcare professional before treatment!



Backstop in the right position during injection

GB INSTRUCTIONS FOR USE FOR BELOTERO® SOFT LIDOCAINE

1. Description

Belotero Soft Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer. Belotero Soft Lidocaine contains 0.3% of lidocaine hydrochloride.

2. Presentation

Belotero Soft Lidocaine is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

3. Composition

Cross-linked sodium hyaluronate: 20.0 mg/ml
Lidocaine hydrochloride: 3.0 mg/ml
Phosphate buffer pH 7 q.s.: 1ml
Sodium hyaluronate is produced by fermentation of *Streptococcus equi*.

4. Indications

4.1. Intended use

Belotero Soft Lidocaine is an injectable biodegradable implant intended for filling of peroral fine lines. The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

4.2. Indications

Belotero Soft Lidocaine is indicated for injection into the superficial to mid-dermis for treatment of perioral fine lines.

5. Posology and administration method

Belotero Soft Lidocaine is designed to be injected into the superficial to mid-dermis by authorized practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection in order to minimize the risk of potential complications.

Belotero Soft Lidocaine can be injected with the blanching technique. Inject Belotero Soft Lidocaine slowly and not too fast to apply the least amount of pressure necessary, according to the appropriate injection technique using the provided needles.

The risk of an intravascular injection can be reduced by different strategies, including aspiration prior to injection, utilizing lower volumes and serial injections in high-risk areas, treating one side at a time, pinching/teasing the skin to provide more space superficial to the branches of the main arteries, and manual occlusion of the origin of the supratrochlear vessels with the non-dominant finger. Blunt cannulas may reduce, but not eliminate the risk.

General recommended injection techniques are for example: linear or serial threading, fanning, cross-hatching or serial puncture. The quantity of product to be injected depends on the area to be corrected. Belotero Soft Lidocaine must be injected into a healthy, non-inflamed skin. Before injection, thoroughly disinfect the area to be treated.

To ensure optimal use of Belotero Soft Lidocaine, it is recommended to assemble the needle according to the diagrams below. Improper assembly may lead to a separation of the needle and syringe and / or leakage of the gel at the Luer lock connection during injection.

If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle.

The quantity of the gel to be injected depends on the area to be treated and the connection to be achieved. Do not over-correct.

The graduations on the syringe label are only intended for orientation for the user.

Gently massage the treated area after the injection to distribute the product uniformly.

Before treatment, the patient's suitability for the treatment and the patient's need for pain relief (topical anesthetics, ice packs, distraction techniques, local anesthetic injections, or nerve blocks) depending on the injection site and size of needle used), should be assessed.

6. Contra-Indications

Belotero Soft Lidocaine is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, lidocaine hydrochloride, BDE or to amide-type local anesthetics;
 - In pregnant and breast-feeding women;
 - In young patients under 18 years old;
 - In patients presenting a general infection;
 - In patients presenting an active auto-immune disease. Do not inject Belotero Soft Lidocaine into blood vessels. Do not inject Belotero Soft Lidocaine into areas presenting cutaneous inflammation or infection due to e.g. immunologic, allergic, bacterial, fungal or viral causes.
- Do not inject Belotero Soft Lidocaine into an area previously treated with a permanent dermal filler.

7. Precautions for use

Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

In the absence of available clinical data on tolerance on the injection of Belotero Soft Lidocaine in patients with anticoagulants or with an active autoimmune disease or in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject Belotero Soft Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

It is recommended not to inject Belotero Soft Lidocaine in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or keloids.

Belotero Soft Lidocaine injected in the glabella, nose, temples and periorbital area may be associated with an increased risk for intravascular complications and the consequences of local vascular occlusion, embolization, vision impairment, blindness, ischemia, necrosis or infarction.

Belotero Soft Lidocaine can be used in combination with other Belotero products during the same session but in different facial areas. Instructions for use of each product should be followed.

No clinical data is available on the injection of Belotero Soft Lidocaine into patient with a Fitzpatrick skin type V/VI.

Belotero Soft Lidocaine can be used in combination treatments such as with botulinum toxin and/or calcium hydroxylapatite (Radiesse®) only if injected in different facial areas. Practitioners should be experienced and patients appropriately selected as benefits but also adverse events can be cumulative and causality of adverse events could become difficult to determine. Instructions for use, depth of injection

and appropriate recommendation of each product should be followed. No clinical data are available on the injection of Belotero Soft Lidocaine into an area already treated with other filling aesthetic products or procedures.

Belotero Soft Lidocaine must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, Belotero Soft Lidocaine must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of Belotero Soft Lidocaine with the above-mentioned treatments.

Patients using anti-coagulation, anti-platelet, or thrombolytic medications (e.g. warfarin), anti-inflammatory drugs (oral/injectable corticosteroids or non-steroidal anti-inflammatory drugs (NSAIDs); e.g. aspirin, ibuprofen), or other substances known to increase coagulation time (vitamins or herbal supplements, e.g. Vitamin E, garlic, Ginkgo biloba and St. John's Wort), from 10 days pre- to 3 days post injection may have increased reactions of hematomas, nodules or bleeding at the injection site.

Injection of Belotero Soft Lidocaine into patients with a history of previous herpetic eruption may be associated with reactivation of the herpes (HHV related diseases, e.g. pityriasis rosea).

In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject Belotero Soft Lidocaine on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.

It should be noted that the presence of lidocaine may cause local redness or hypersensitivity or transient loco-regional numbness.

For normal healthy adults, it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in signs of the central nervous system or cardiovascular toxicity.

When used concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or agents structurally related to amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.

Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents.

Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the inner packaging has been opened or damaged.

Do not transfer Belotero Soft Lidocaine into another container and do not add other substances to the product.

Only the gel is sterile, but not the outside of the syringe. Discard the syringe and the remaining product and the needles in the appropriate container after use.

Do not re-sterilize and do not reuse due to the associated risks including infection.

The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area and should avoid strenuous physical activity following treatment.

The patient must avoid drinking alcohol for 24 hours before and after treatment. Alcohol may cause the blood vessels to dilate and cause more bruising.

paralysis, facial paralysis, trigeminal neuralgia).

Literature reported rare cases of the following adverse events with hyaluronic acid products such as infection (e.g. erysipelas, phlegmon, cellulitis, including open or draining wounds and (dental) abscess, impetigo, pustules) chronic infection (including biofilm formation), scarring, persistent skin discoloration, sensory dysfunction, infarction, non-thrombotic lung embolism as well as sarcoïd granuloma formation in subjects with hepatitis C and interferon treatment, cerebral injuries (e.g. intracranial penetration, subarachnoid hemorrhage), strabismus, ophthalmoplegia, iris adhesions, cataract, conjunctival hemorrhage, eyelid ptosis and lacrimation.

The risk of granuloma, ischemia, necrosis and vascular occlusion is higher with deep injections and high volumes. Isolated cases of visual impairment or blindness following unintentional intra-arterial injection have been reported in literature.

Patients should be instructed to report any side effects which last for more than one week and any adverse event as soon as it occurs to his/her practitioner, especially if patient has changes in his/her vision, signs of a stroke (including sudden difficulty speaking, numbness or weakness in his/her face, arms or legs, difficulty walking, face drooping, severe headache, dizziness, or confusion), white appearance of the skin, or unusual pain during or shortly after treatment. The practitioner may then refer the patient to the appropriate treatment.

Patients with lighter skin types are more likely to develop injection-related adverse events. However, patients with skin of color are more likely to develop post-inflammatory hyperpigmentation and / or hypertrophic scar/keloid formation following injection procedures. Patients with specific ethnic characteristics, e.g. Asian population, should be informed of a higher risk of tissue reactions, e.g. itching, swelling, erythema, inflammation.

10. Assembly of needle to syringe
For optimal use of Belotero Soft Lidocaine, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

- (1) Firmly hold the glass cylinder of the syringe and the Luer-lock adaptor between the thumb and the index fingers.
- (2) Grasp the protective cap with the other hand and unscrew it.
- (3) Push & Twist the needle on the syringe until a resistance is felt. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
- (4) Keep holding the Luer-lock and remove the sheath from the needle.

11. Storage
Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

12. References
Updated documentation may be available from ANTEIS SA, Switzerland.